



Quick Protocol

Ovation Kit
Part Number

Ovation Kit
Lot Number

Date

Generation of ST-cDNA (sense-target cDNA)

Obtain all **Exon Module Reagents** from -20°C.

Invert **E3** vial to mix and spin. Thaw **E1** and **E2**, vortex, spin, and place on ice.

Pipette 3 µg SPIA™ cDNA into a 0.2 mL PCR tube for each reaction. Bring volume up to 20 µl with water if necessary.

Add 6 µl **E1** Primer Mix to each sample and mix well.

Place the tubes in a thermal cycler running program 1 (95°C - 5 min, 4°C - forever).

Make **ST-cDNA Master Mix**. Per sample combine:

38 µl **Buffer Mix E2** and 6 µl **Enzyme Mix E3**.

Enter the number of samples and volumes of each reagent into this table.

No. of rxns	E2	E3
1	38 µl	6 µl

Mix the **ST-cDNA Master Mix**, spin and place on ice.

For each reaction, add 44 µl **ST-cDNA Master Mix** and mix well.

Place the tubes in a thermal cycler running program 2 (4°C - 1 min, 30°C - 10 min, 42°C - 60 min, 75°C - 10 min, 4°C - forever).

Once the thermal cycler reaches 4°C, remove tubes and spin.

Proceed to purification, fragmentation and labeling or store ST-cDNA at -20°C.

Notes

